



**Pariksha
Biotech**
A game changer in IVD

ALP Estimation Kit

High-Q Alkaline Phosphatase-ML(Ultra)

IFCC-AMP Buffer Method




Intended Use:

In vitro test for the quantitative determination of alkaline phosphatase (ALP) in human serum.

Clinical Significance

Alkaline phosphatase in serum consists of four structural genotypes: the liver-bone-kidney type, the intestinal type, the placental type and the variant from germ cells. It occurs in osteoblasts, hepatocytes, kidneys, spleen, placenta, prostate, leukocytes and the small intestine. The liverbone-kidney type is particularly important. A rise in the Alkaline Phosphatase activity occurs with all forms of cholestasis, particularly with obstructive jaundice. It is also elevated in diseases of the skeletal system, such as Paget's disease, hyperparathyroidism, rickets and osteomalacia, as well as with fractures and malignant tumors. A considerable rise in the Alkaline Phosphatase activity is sometimes seen in children and juveniles. It is caused by increased osteoblast activity following accelerated bone growth. Various reference values for the purposes of clinical evaluation have been assigned to differing age groups. In 1946, Bessey, Lowry and Brock published a method for the determination of Alkaline Phosphatase using p-nitrophenyl phosphate as substrate buffered with glycine/NaOH. This assay described here meets the IFCC recommendations and uses AMP buffer.

PRINCIPLE :



Active Ingredients of reagents:

AMP	pH 9.7	3.5 Mol/l
Magnesium chloride		0.6 mMol/l
p-Nitrophenylphosphate Analogue		0.070 mMol/l
Sodium Azide		0.10 %
Zinc Sulphate		0.3 mMol/l

Reagent Storage and Stability

- When the reagent is stored properly at 2-8°C and the contamination is avoided, it is stable up to the expiry date mentioned on the labels and kit box.
- It is recommended that when the reagent is not in use, it should be strictly kept at 2-8°C for better performance and stability.
- Substrate Reagent must be protected from direct sun light.
- Do not expose the reagent to high temperatures above 10°C.

Specimen

- Unhemolysed serum is the preferred specimen. DO NOT USE PLASMA
- Serum samples are stable for 24 hours at 2-8°C and for a month at -20°C. Please note that ALP activity in serum stored at 2-8°C increases with time.

Test Procedure:

High-Q Alkaline Phosphatase is a Mono Liquid Reagent and is ready to use. There is no reconstitution required.

Test Procedure:

Substrate Reagent	500 µl
Serum	10 µl

Mix well and aspirate in to the analyzer. After 20 Sec incubation (Delay) measure the change of optical density during the next 90 seconds (Measuring Time) against distilled water at 405 nm as follows.

Ao-	Exactly after 20 Seconds.
A1	Exactly after 90 seconds

Calculation

Calculate the average change in absorbance per minute ($\Delta A/\text{min}$) and multiply by the corresponding factor.

ALP activity [IU/L]: $\Delta A/\text{min} \times 3500$

System Parameters:

Reaction Type (Mode)	:	Kinetic
Reaction Direction	:	Increasing
Wave Length	:	405 nm
Flow Cell Temp.	:	37°C
Zero Setting With	:	Distilled Water
Delay Time	:	20 Seconds
Measuring Time	:	90 Seconds
Reagent Volume	:	500 µl
Sample Volume	:	10 µl
Factor	:	3500
Linearity	:	1250
Units	:	IU/L

Expected Values:

As per the recommendations of International Liver Foundation the following reference values are assigned for ALP IFCC Assays. But it is for the Laboratory to establish its own reference values.

Reference Range:

Adults :

Men	[IU/L]	40 - 130
Women	[IU/L]	35 - 105

Children:

	Female	Male
1-30 Days	[IU/L] 48-406	75-316
1 Month- 1 Year	[IU/L] 124-341	82-383
1 Year- 3 Years	[IU/L] 108-317	104-345
4 Years-6Years	[IU/L] 96-297	93-309
7 Years-9 Years	[IU/L] 69-325	86-315
10 Years- 12 Years	[IU/L] 51-332	42-362
13 Year- 15 Years	[IU/L] 50-162	74-390
16 Year- 18 Years	[IU/L] 47-119	52-171

Specificity /Interference

No interference is absorbed by triglycerides up to 2200 mg/dl, bilirubin up to 42 mg/dl, hemoglobin up to 150 mg/dl, intra lipids up to 1200 mg/dl and ascorbic acid up to 32 mg/dl

Sensitivity

The lower limit of detection is 2 IU/L



**Pariksha
Biotech**
A game changer in IVD

ALP Estimation Kit



High-Q Alkaline Phosphatase-ML(Ultra)

IFCC-AMP Buffer Method

Performance Characteristics:

Linearity - Measuring Range

1250 IU/L

The test has been developed to determine Alkaline Phosphatase activities which correspond to maximal A/min of 0.227

If such value is exceeded the sample should be diluted 1+9 with NaCl solution (9g/l) and the results be multiplied by 10.

Precision

Intra-assay precision

N=20	mean (IU/L)	SD (IU/L)	CV (%)
Sample 1	143.1	2.30	1.59
Sample 2	448.4	6.60	1.46
Sample 3	449.2	5.90	1.32

Method Comparison

A comparison of the High-Q ALP - ML (y) with a commercial obtainable assay (x) gave following results with 80 samples.

$$y = 0.992x + 2.64; r = 0.991$$

Notes :

- Avoid contact with skin and mucous membranes. The reagents contain sodium azide as preservative. Do not swallow.
- During the reaction p-nitrophenol is produced. This substance is poisonous when inhaled, swallowed or when absorbed through the skin. If the reaction mixture comes into contact with skin or membranes, wash copiously with water and consult a doctor.

References

- Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.
- Bessey OAH et al. J Biol Chem 1946;164:321
- Empfehlungen der Deutschen Gesellschaft für Klinische Chemie. Standard-Methode zur Bestimmung der Aktivität der alkalischen Phosphatase. Z klin Chem u klin Biochem 1972;10:191.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-474.

Ordering information

Ref./Cat.

Pack Size

Presentation

P-ALP(ML) - 50	5 x 10 ml	Single Liquid Reagent
P-ALP(ML) - 100	10 x 10 ml	

Product Features

- Liquid Stable, Ready to use Mono Reagent
- PNPP Analogue used for better stability
- 110 Seconds increasing kinetic reaction (20 Sec Delay+ 90 Sec Measuring)
- Measuring Wavelength 405 nm
- Kinetic factor 3500 at 37° C.
- Linearity 1250 IU/L
- Serum is the only specimen
- Available as multipurpose reagents and dedicated system packs



Symbols used with IVD devices

	Date of manufacture		Manufactured by
	In vitro diagnostic medical device		Keep away from sunlight
	Do not freeze		This way up
	Use-by date		Do not use if package is damaged
	Keep dry		Batch Code
	Consult instructions for use		Reagent
	Consult electronic instructions for use		Calibrator Material
	Catalogue number		
	Caution		

Manufactured by :-
Pariksha Biotech Pvt Ltd,
Plot no.1/B-14,SVICE,
Balanagar,
Hyderabad-500037
Telangana State
CIN No: U24232TG2011PTC072118



Pariksha's world inside
SCAN TO EXPLORE MORE



+91 7075706709



info@parikshabio.com



www.parikshabio.com