



**Pariksha
Biotech**
A game changer in IVD

High-Q Procalcitonin (PCT)

Latex Enhanced Turbidimetric Immuno Assay (LETIA)

Intended Use: Kit for the In Vitro quantitative determination of Procalcitonin in Human Serum or Plasma



Intended Use

High-Q Procalcitonin Reagent Kit is a Latex-Enhanced Turbidimetric Assay intended for *in vitro* quantitative detection of Procalcitonin in human serum or plasma on automated clinical chemistry analyzers.

Summary and Explanation

Procalcitonin (PCT) is a 116 amino acid protein, the prohormone of calcitonin. Whereas hormonally active calcitonin is produced exclusively in the C-cells of the thyroid gland after specific intracellular proteolytic processing of the prohormone PCT, PCT is ubiquitously and uniformly expressed in multiple tissues throughout the body in response to sepsis¹. In healthy conditions, the PCT levels in the circulation are very low (< 0.05 ng/ml). Elevated circulating levels of PCT are important indicators in response to microbial infections and a powerful tool in the early detection of sepsis. According to the criteria of the American College of Chest Physicians/Society of Critical Care Medicine Consensus Conference, rediagnosis of systemic bacterial infection sepsis is categorized as follows: PCT < 0.5 ng/mL: systemic infection (sepsis) is not likely, local bacterial infection is possible; PCT ≥ 0.5 ng/mL and < 2 ng/mL: systemic infection (sepsis) is possible but other conditions known to increase PCT as well; PCT ≥ 2 ng/mL and < 10 ng/mL: systemic infection (sepsis) is likely, unless other cause are known; PCT ≥ 10 ng/mL, important systemic inflammatory response, almost exclusively due to severe bacterial or septic shock. Deferential diagnosis of lower respiratory tract infections: PCT < 0.25 ng/mL: bacterial infection unlikely. Use of antibiotic is discouraged; PCT ≥ 0.25 ng/mL and < 0.5 ng/mL: bacterial infection is possible. Recommendation to initiate antimicrobial therapy; PCT ≥ 0.5 ng/mL: suggests the presence of bacterial infection. Antibiotic treatment strongly recommended.

Test Principle

The Kit utilizes latex-enhanced immunoturbidimetry to measure the PCT level in human serum or plasma. During the test, PCT in the sample binds with the antibody that is coated on latex particles to cause agglutination. The turbidity caused by agglutination is detected optically by chemistry analyzer at wavelength of 630 nm. The change in absorbance is proportional to the level of PCT in the sample. The actual concentration is obtained by comparing with a calibration curve with known concentrations.

Reagents:

R1	Tris Buffer solution	(0.2 mol/L)
	NaCl	(6.0 g/L)
	PEG-6000	(22 g/L)
R2	Suspension of anti-human PCT monoclonal antibody coated latex particles	(0.3%)
	Phosphate Buffered Saline (ph=7.00)	(0.2 mol/L)
	BSA	(3 g/L)
	Tween20	(3 mL/L)

- PCT calibrator set is provided along with the Reagent Kit PCT Control set is purchased optional.

Specimen Collection and Handling

Follow standard laboratory procedures to collect serum or lithium heparin and EDTA plasma samples. It is recommended to perform test immediately after sample collection to avoid hemolysis. If the test cannot be done immediately, store sample at 4° C for up to 3 days or at -20° C for up to 6 months. Avoid repeated freezing and thawing.

Usage and Instrumentation

The Kit is applicable on most automated chemistry analyzers. Refer to specific instrument application sheets

Storage and Stability:

Store the reagents at 2-8°C. Avoid direct sunlight. The Kit is stable till the expiry date mentioned on the labels when stored properly

Precautions

- The Kit is for *in vitro* diagnostic use only.
- The procedural instructions must be followed to obtain accurate results.
- Do not use the reagents beyond the expiration date.
- Do not mix and use different lots of reagents.
- Treat all specimens as infectious. Proper handling and disposal procedures of specimens and test materials should be strictly followed.
- Reagents contain sodium azide as preservative. Avoid contact with skin and eyes, flush with copious amounts of water when disposing.

Limitations:

The Kit is for *in vitro* use on automated chemistry analyzers only. Hemolysis samples may cause inconsistent results. The test result from the Kit should not be used as the only basis for definite diagnosis.

Quality Control:

High-Q PCT control is recommended before each batch of tests to ensure the test is properly performed and all reagents and the instrument are functional as specified.



**Pariksha
Biotech**
A game changer in IVD



**Pariksha
Biotech**
A game changer in IVD

High-Q Procalcitonin (PCT)

Latex Enhanced Turbidimetric Immuno Assay (LETIA)

Intended Use: Kit for the In Vitro quantitative determination of Procalcitonin in Human Serum or Plasma



Semi Auto Analyzers Assay Procedure: (Multi Point Calibration with 5 Calibrator Levels)

System Parameters:

Reaction Type (Mode)	Fixed Time- Non Linear- Multi Standard-Spline
Reaction Direction	Increasing
Wave Length	630 nm (600-630 nms)
Flow Cell Temp.	37°C
Delay Time	60 Seconds
Measuring Time	120 Seconds
Blank	Distilled Water Blank
Reagent Volume	360 µl (R1) + 120 µl (R2)
Sample Volume	50 µl
Calibrator Concentrations	(On the Vials Lot Specific)
Linearity	60 ng/mL

Procedure :

Reagent	Calibrator	Serum/Control
R1	360 µl	360 µl
Calibrator	50 µl	----
Serum/Control	----	50 µl
Mix and incubate for 5 Minutes at 37 °C		
R2	120 µl	120 µl

- 1) Read absorbance A1 after 60 Seconds. (Delay)
- 2) Incubate and Read the absorbance A2 after 120 Seconds (Measuring)
- 3) Calculate the absorbance differences $\Delta A = A2 - A1$ for each point of the calibration curve, controls and all unknown samples.
- 4) The concentration of PCT in the unknown sample can be calculated from $\Delta A = A2 - A1$
- 5) Using a 3rd order polynomial mathematical model where abscissa (X) is the $\Delta A = A2 - A1$ and ordinate (Y) is the concentration of PCT or plotting the values of $\Delta A = A2 - A1$ obtained for every concentration level of the calibrator against the PCT concentration and interpolating the individual $\Delta A = A2 - A1$ of every sample in the calibration curve.

Calculations with Calibrators/ Calibration Curve/ Result Interpretation:

Calculation

The concentration of PCT in unknown samples is derived from a calibration curve using an appropriate mathematical model such as spline. The calibration curve is obtained with 5 calibrators at deferent levels.

Stability of calibration: 4 weeks

Performance Characteristics:

Linearity: 0.05-60 ng/mL ($R^2 \geq 0.990$)

Accuracy: Relative deviation $\leq 15\%$

Precision: Within Run: $CV \leq 10\%$
Run-to-Run: $CV \leq 10\%$

Interference: no interference detected for: Triglyceride (1000 mg/dL), Ascorbic Acid (10mM), Hemoglobin (400 mg/dL), billirubin (40 mg/dL), and RF (200 IU/mL)

Expected Values:

Less than 0.5 ng/mL

It is recommended that each laboratory should establish it's own expected range

Result Interpretation:

According to the criteria of the American College of Chest Physicians/Society of Critical Care Medicine Consensus Conference, diagnosis of systemic bacterial infection sepsis is categorized as follows:

PCT < 0.5 ng/mL:

Systemic infection (Sepsis) is not likely, Local bacterial infection may possible.

≥ 0.5 ng/mL and ≤ 2 ng/mL:

PCT ≥ 0.5 ng/mL and ≤ 2 ng/mL: Systemic infection is possible, but other non bacterial conditions are known to increase PCT as well. Should be clinically correlated before starting the Antibiotic Treatment

PCT ≥ 2 ng/mL and < 10 ng/mL:

Systemic infection (sepsis) is likely, unless other causes are known; High risk for progression to severe systemic infection (severe sepsis).

PCT ≥ 10 ng/mL,

important systemic inflammatory response, due to severe bacterial or septic shock.

Clinical Limitations and Correlations

Increased PCT levels may not always be related to systemic bacterial infection

Several situations have been described where PCT can be elevated by non-

bacterial causes. These include, but are not limited to:

- neonates < 48 hours of life (physiological elevation)
- the first days after a major trauma, major surgical intervention severe burns, treatment with OKT3 antibodies and other drugs stimulate the release of pro-inflammatory cytokines

Patients with invasive fungal infections, acute attacks of Plasmodium falciparum malaria

Patients with prolonged or severe cardiogenic shock, prolonged severe organ perfusion anomalies, small cell lung cancer, medullary C-cell carcinoma of the thyroid.

Low PCT levels do not automatically exclude the presence of bacterial infection.

Such low levels may be obtained, during the early course of infections, in localized infections and in subacute endocarditis. Therefore, follow-up and re-evaluation of PCT in clinical suspicion of infection is pivotal. The PCT measuring technique should be chosen according to clinical use.

Final clinical diagnosis based on the test result must be supported by a comprehensive judgment of the concerned physician including clinical symptoms and other relevant test results.

Intended Use: Kit for the In Vitro quantitative determination of Procalcitonin in Human Serum or Plasma

Fast and highly specific PCT increase in bacterial infection and sepsis

One major advantage of PCT compared to other parameters is its early and highly specific increase in response to severe systemic bacterial infections and sepsis. Therefore, in septic condition increased PCT levels can be observed 3-6 hours after an infectious challenge.

PCT levels are usually low in viral infections chronic inflammatory disorders or autoimmune disorders. PCT levels in sepsis are generally greater than 0.5-2 ng/mL and often reach values between 10 and 10 ng/mL, or considerably higher in individual cases, thereby enabling diagnostic differentiation between these various clinical conditions and a severe bacterial infection (sepsis)

DEFINITIONS

Definitions for the terms of "SIRS", "sepsis", "severe sepsis" or "septic shock" have been proposed by the ACCP/SCCM consensus conference in 1992, and are now widely used (see below table 1).

SIRS and sepsis definition (ACCP/SCCM-criteria) SIRS (Systemic Inflammatory Response Syndrome)	2 or more of the following criteria: • Temperature > 38°C or < 36°C • Heart rate > 90 beats/min • Respiratory rate > 20 breaths/min or PaCO ₂ < 32 mm Hg (<4.3 kPa) • WBC > 12 000 cells/μL or < 4 000 cells/μL or > 10% immature (band) forms
Sepsis	Documented infection together with 2 or more SIRS criteria
Severe Sepsis	Sepsis associated with organ dysfunction including, but not limited to, lactic acidosis, oliguria, hypoxemia, coagulation disorders, or an acute alteration in mental status.
Septic Shock	Sepsis with hypotension, despite adequate fluid resuscitation along with the presence of perfusion abnormalities patients who are on inotropic or vasopressor agents may not be hypotensive at the time when perfusion abnormalities are detected.

ACCP: American College of Chest Physicians
SCCM: Society of Critical Care Medicine

[REFERENCES]
















- Müller B, *et al.*, Ubiquitous expression of the calcitonin-receptor-like receptor 1 gene in multiple tissues in response to sepsis. *J Clin Endocrinol Metab* 2001; 86(1):396-404.
- Meisner M. Procalcitonin (PCT) – A new, innovative infectious parameter. Biochemical and clinical aspects. Thieme; Stuttgart, New York, 2000; ISBN 3-13-105503-0.
- Christ-Crain M, *et al.*, Procalcitonin in bacterial infection – hype, hope or more or less? *Swiss Med Wkly* 2005; 135: 451-60.

High-Q Procalcitonin

Product Features

- ❖ Latex Enhanced Turbidimetric Immuno Assay (LETIA)
- ❖ Liquid Stable Two Reagents
- ❖ 5 Level Calibrator Set provided
- ❖ Measurement at 630 nms (600-630 nms)
- ❖ Test Procedure time 8 minutes at 37°C
- ❖ Linearity : 60.0 ng/mL
- ❖ Adaptable to Semi and Automated Analyzers

Symbols used with IVD devices

	Date of manufacture		Manufactured by
	In vitro diagnostic device		Keep away from sunlight
	Do not freeze		This way up
	Use by (yyyy-mm-dd or mm/yyyy)		Reagent
	Calibrator Material		Batch code
	Temperature limitation (store at)		Control
	Consult instructions for use		Keep dry Keep away from rain
	Catalog Number		

